

## PHYSICAL AND CHEMICAL PROPERTIES

**Product Name:** Ganoderic acid Epsilon

**Catalog Number:** BCX1480

**Cas No.:** 294674-05-8

**Formula:** C<sub>30</sub>H<sub>44</sub>O<sub>7</sub>

**M.Wt:** 516.67

**Purity:** >98%


**Physical Description:** N/A

**Synonyms:** Ganoderic acid ε

**Chemical Name:** (E,4S,6R)-6-[(3S,5R,7S,10S,13R,14R,17R)-3,7-dihydroxy-4,4,10,13,14-pentamethyl-11,15-dioxo-2,3,5,6,7,12,16,17-octahydro-1H-cyclopenta[a]phenanthren-17-yl]-4-hydroxy-2-methylhept-2-enoic acid

**SMILES:** CC(CCC(C=C(C)C(=O)O)O)C1CC(=O)C2(C1(CC(=O)C3=C2C(CC4C3(CCC(C4(C)C)O)C)O)C)C

**Solubility:** N/A



## STORAGE AND SHIPPING

**Storage:** Store the product in sealed, cool and dry condition.

**General tips:** For obtaining a higher solubility, please warm the tube at 37° and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20° for several months.

**Shipping Condition:** Ship with RT , or blue ice upon request

## INTENDED USE

1. Reference standards;
2. Pharmacological research;
3. Food research;
4. Synthetic precursor compounds;
5. Intermediates & Fine Chemicals;
6. Others.

## CAUTION!

**FOR RESEARCH PURPOSES ONLY.**

**NOT FOR HUMAN, VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.**

*Specific storage and handling information for each product is indicated on the product datasheet. Most BioCrick products are stable under the recommended conditions. Products are sometimes shipped at a temperature that differs from the recommended storage temperature. Short-term storage of many products are stable in the short-term at temperatures that differ from that required for long-term storage. We ensure that the product is shipped under conditions that will maintain the quality of the reagents. Upon receipt of the*

*product, follow the storage recommendations on the product data sheet.*

---

F9-15, Building B2, Tianfu Life Science Park, 88 South Keyuan Road, Hi-Tech Zone, Chengdu, Sichuan 610041, PRC  
Tech Support: service@biocrick.com Tel: +86-28-8543-3893 FAX: +86-28-8543-3893